

Human S100A11 Recombinant Protein



RPPB4548

Product Information Protein Information

Product SKU:

RPPB4548

Accession:

P31949

Host:

Escherichia Coli.

Protein description:

S100A11 Human Recombinant is expressed in E. coli having a molecular weight of 17kDa fused to an amino terminal hexahistidine tag.

Synonyms:

Protein S100-A11, S100 calcium-binding protein A11, Calgizzarin, MLN 70, S100A11, MLN70, S100C.

Formulation:

S100A11 is supplied in PBS and 50% glycerol.

Purity:

Greater than 95% as determined by SDS-PAGE. 2 bands on Western blot at 17 and 34 kDa, respectively representing monomeric and dimeric form.

Stability:

Store at 4°C if entire vial will be used within 1-2 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.