

Protein A Recombinant Protein



RPPB4787

Product Information Protein Information

Product SKU:

RPPB4787

Accession:

P38507

Host:

Escherichia Coli.

Protein description:

Recombinant Staphylococcal Protein A produced in E.Coli is a non-glycosylated, Polypeptide chain having a molecular mass of 45 kDa. Recombinant Staphylococcal Protein A is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear colorless solution.

Synonyms:

Immunoglobulin G-binding protein A, IgG-binding protein A, Staphylococcal protein A, SPA.

Formulation:

The protein solution contains no additives.

Purity:

Greater than 95.0% as determined by RP-HPLC.

Stability:

SPA should be stored at -20°C.

Biological Activity:

Greater than 95.0% binding activity to human IgG.