

# T.Vaginalis P65 Recombinant Protein



RPPB4893

---

## Product Information    Protein Information

**Product SKU:**

RPPB4893

**Host:**

Escherichia Coli.

**Protein description:**

The Recombinant Trichomonas vaginalis p65, E.Coli derived migrates around 67kda on SDS-PAGE. The full length protein is purified by standard chromatography techniques.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

T. Vaginalis p65protein solution in PBS and 2M Urea.

**Purity:**

Protein is >90% pure as determined by 12% PAGE (coomassie staining).

**Stability:**

T. Vaginalis p65 Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.