

Human TAF15 Recombinant Protein



RPPB4900

Product Information Protein Information

Product SKU:

RPPB4900

Accession:

Q92804

Host:

Escherichia Coli.

Protein description:

TAF15 Human Recombinant produced in E. coli is a single polypeptide chain containing 282 amino acids (148-406) and having a molecular mass of 30kDa. TAF15 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Npl3, RBP56, TAF2N, TAFII68, TATA-binding protein-associated factor 2N, 68 kDa TATA-binding protein-associated factor, RNA-binding protein 56, TAF15.

Formulation:

The TAF15 solution (0.25mg/ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 30% glycerol and 1mM DTT.

Purity:

Greater than 90% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MGSSHHHHHH SGLVPRGSH MGSSYHSQRE NYSHHTQDDR RDVSRYGEDN RGYGGSQGGG
RGRGGYDKDG RGPMTGSSGG DRGGFKNFGG HRDYGPRTDA DSES DNSDNN TIFVQGLGEG VSTDQVGEFF
KQIGIIKTNK KTGKPMINLY TDKDTGKPKG EATVSFDDPP SAKAAIDWFD GKEFHGNIK VSFATRRPEF
MRGGGSGGGR RGRGGYRGRG GFQGRGGDPK SGDWWCPNPS CGNMNFARRN SCNQCNEPRP
EDSRPSGGDF RGRGYGGERG YR.