

# Human TAX1BP3 Recombinant Protein



RPPB4911

---

## Product Information    Protein Information

**Product SKU:**

RPPB4911

**Accession:**

O14907

**Host:**

E.coli.

**Protein description:**

TAX1BP3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 144 amino acids (1-124 a.a.) and having a molecular mass of 15.8kDa.TAX1BP3 is fused to a 20 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered colorless solution.

**Synonyms:**

Tax1 (human T-cell leukemia virus type I) binding protein 3, TIP-1, Glutaminase-interacting protein 3, Tax interaction protein 1.

**Formulation:**

TAX1BP3 protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH8.0), 1mM DTT, 0.1M NaCl and 20% glycerol.

**Purity:**

Greater than 95% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles.

**Amino Acid Sequence:**

MGSSHHHHHH SSGLVPRGSH MSYIPGQPVT AVVQRVEIHK LRQGENLILG FSIGGGIDQD PSQNPFSEDK  
TDKGIYVTRV SEGGPAEIAG LQIGDKIMQV NGWDMTMVTH DQARKRLTKR SEEVRLLVLT RQSLQKAVQQ  
SMLS.