

Bovine Sm Recombinant Protein



RPPB4941

Product Information Protein Information

Product SKU:

RPPB4941

Host:

Bovine tissues.

Protein description:

Bovine Small Nuclear Ribonucleoprotein Polypeptides purified from bovine tissues by proprietary protein-chemical techniques.

Appearance:

Sterile Filtered clear solution.

Formulation:

SNRPD is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.

Purity:

Greater than 60% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.