

Bovine Trypsin Recombinant Protein



RPPB5025

Product Information Protein Information

Product SKU:

RPPB5025

Host:

Corn.

Protein description:

Recombinant Bovine Trypsin is free from any animal and human sources. Trypsin Bovine specifically cleaves peptide bonds after basic amino acids such as lysine and arginine.

Appearance:

Sterile Filtered lyophilized powder.

Formulation:

The protein was lyophilized without any additives.

Purity:

Greater than 90% as determined by SDS-PAGE.

Solubility:

It is recommended to reconstitute the lyophilized Bovine Trypsin in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Stability:

Store the Bovine Trypsin between 2-8°C, do not freeze.

Biological Activity:

4,313 Units/mg.