

U. Urealyticum Recombinant Protein



RPPB5061

Product Information Protein Information

Product SKU:

RPPB5061

Host:

E.Coli

Protein description:

MB antigen of Ureaplasma parvum serovar fused with GST, is expressed in E.Coli, having a Mw of 40kDa, and recognized by the specific IgM and IgG to Ureaplasma parvum.

Appearance:

Sterile filtered colorless solution.

Synonyms:

U. urealyticum, UreD, MB Antigen.

Formulation:

MB Antigen (2.2 mg/ml) is formulated in 1x PBS and 1mM EDTA pH 7.4 and 0.01% Sodium azide.

Purity:

Greater than 95% as determined by 10% PAGE (Coomassie staining).

Stability:

Upon arrival, Store at -20°C. Please prevent freeze-thaw cycles.