

Human RPS27A Recombinant Protein



RPPB5063

Product Information Protein Information

Product SKU:

RPPB5063

Accession:

P62988

Host:

Escherichia Coli.

Protein description:

Ubiquitin Human Recombinant expressed in E.coli and purified by ion-exchange chromatography, contains 76 amino acids and has an Mw of 8.6kDa.

Appearance:

Sterile Filtered solution.

Synonyms:

Ubiquitin, Ribosomal Protein S27a, CEP80, UBA80, UBCEP1, UBCEP80, HUBCEP80, RPS27A.

Formulation:

1mg/ml containing 50mM HEPES (pH7.5) 150mM NaCl and 10% glycerol.

Purity:

Greater than 95% as estimated by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MQIFVKLTG KTITLVEPS DTIENVKAKI QDKEGIPPDQ QRLIFAGKQL EDGRTLSDYN IQKESTLHLV LRLRGG.