

Human Ubiquitin K48R Recombinant Protein



RPPB5066

Product Information Protein Information

Product SKU:

RPPB5066

Accession:

P62988

Host:

Escherichia Coli.

Protein description:

Recombinant human ubiquitin featuring a Lys 48 to Arg48 mutation is useful for the reduction of poly-Ub chain length and conjugation rates. Ubiquitin K48R is expressed in E.coli and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Ubiquitin, Ribosomal Protein S27a, CEP80, UBA80, UBCEP1, UBCEP80, HUBCEP80, RPS27A, Ubiquitin K48R.

Formulation:

Diluted in PBS plus 5% glycerol.

Purity:

Greater than 95% as determined by SDS-PAGE.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.