

Human UBTD2 Recombinant Protein



RPPB5072

Product Information Protein Information

Product SKU:

RPPB5072

Accession:

Q8WUN7

Host:

E.coli.

Protein description:

UBTD2 Human Recombinant produced in E. coli is a single polypeptide chain containing 257 amino acids (1-234) and having a molecular mass of 28.6 kDa. UBTD2 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered colorless solution.

Synonyms:

Ubiquitin domain containing protein 2, Dendritic cell-derived ubiquitin-like protein, Ubiquitin-like protein SB72, DCUBP, MGC30022.

Formulation:

The UBTD2 solution (0.5mg/1ml) contains 20mM Tris-HCl buffer (pH 8.0), 0.1M NaCl, 1mM DTT and 20% glycerol.

Purity:

Greater than 95% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MGSSHHHHHH SSGLVPRGSH MGSMGGCVGA QHDSSGSLNE NSEGTGVALG RNQPLKKEKP
KWKSDYPMTD GQLRSKRDEF WDTAPAFEGR KEIWDALKAA AHAFESNDHE LAQAIIDGAN ITLPHGALTE
CYDELGNRYQ LPVYCLAPPI NMIEEKSDIE TLDIPEPPPN SGYECQLRLR LSTGKDLKLV VRSTDVTFHM
KRR LHAAEGV EPGSQRWFFS GRPLTDKMKF EELKIPKDYV VQVIVSQPVQ NPTPVEN