

Human UTP23 Recombinant Protein



RPPB5098

Product Information Protein Information

Product SKU:

RPPB5098

Accession:

Q9BRU9

Host:

E.coli.

Protein description:

UTP23 Human Recombinant produced in E. coli is a single polypeptide chain containing 272 amino acids (1-249) and having a molecular mass of 30.8 kDa. UTP23 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered colorless solution.

Synonyms:

UTP23 Small Subunit (SSU) Processome Component, Homolog (Yeast), Chromosome 8 Open Reading Frame 53, RRNA-Processing Protein UTP23 Homolog, MGC14595, C8orf53.

Formulation:

The UTP23 solution contains 20mM Tris-HCl buffer (pH 8.0), 0.2M NaCl, 1mM DTT and 20% glycerol.

Purity:

Greater than 90% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MGSSHHHHHH SGLVPRGSH MGSMKTRQK HAKKHLGFFR NNFGVREPYQ ILLDGTFCQA ALRGRIQLRE
QLPRYLMGET QLCTTRCVLK ELETGKDLY GAKLIAQKCQ VRNCPHFKN VSGSECLLSM VEEGNPHHYF
VATQDQNLVS KVKKKPGVPL MFIIQNTMVL DKPSPKTIAF VKAVESGQLV SVHEKESIKH LKEEQGLVKN
TEQSRKRRK KISGPNPLSC LKKKKKAPDT QSSASEKKRK RKRIRNRSNP KVLSEKQNAE GE