

Thrombin Recombinant Protein



RPPB5143

Product Information Protein Information

Product SKU:

RPPB5143

Protein description:

Thrombin was purified from human plasma.

Host:

Human Plasma.

Appearance:

Sterile filtered liquid formulation.

Formulation:

The Factor-IIa contains 50mM sodium Citrate, 0.2M NaCl & 0.1% PEG-8000 pH 6.5.

Stability:

Store at 4°C if entire vial will be used within 7 days. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Biological Activity:

The specific activity was found to be 3120 NIH Units/mg.