

# Human Thyroglobulin Recombinant Protein



RPPB5150

---

## Product Information   Protein Information

**Product SKU:**

RPPB5150

**Accession:**

P01266

**Host:**

Native, Isolated from  
human thyroid  
glands.

**Protein description:**

Human Thyroglobulin is a glycosylated, polypeptide chain having a total molecular mass of 662 kDa (331 kDa per subunit).

**Appearance:**

Sterile Filtered solution.

**Synonyms:**

Thyroglobulin, TGN, AITD3, TG.

**Formulation:**

Human Thyroglobulin is supplied in 20mM HEPES buffer pH-7.6, 150mM NaCl and 40% Sucrose (w/v).

**Purity:**

Greater than 95% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.