

Human Thyroglobulin Recombinant Protein



RPPB5151

Product Information Protein Information

Product SKU:

RPPB5151

Accession:

P01266

Host:

Native, Isolated from human thyroid glands.

Protein description:

Human Thyroglobulin is a biotinylated, glycosylated, polypeptide chain having a total molecular mass of 662 kDa (331 kDa per subunit).

Appearance:

Sterile Filtered solution.

Synonyms:

Thyroglobulin, TGN, AITD3, TG.

Formulation:

Human Thyroglobulin biotinylated is supplied at a 20mM HEPES buffer pH-7.6, 150mM NaCl and 40% Sucrose (w/v).

Purity:

Greater than 80% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.