

# Y. Enterocolitica (O:9) LcrV Recombinant Protein



RPPB5166

---

## Product Information    Protein Information

**Product SKU:**

RPPB5166

**Host:**

Escherichia Coli.

**Protein description:**

Recombinant Yersinia Enterocolitica (O:9) LcrV produced in E.coli is a non-glycosylated, polypeptide chain having a calculated molecular mass of 38,668 Dalton. Y. Enterocolitica (O:9) LcrV is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Formulation:**

Y. Enterocolitica (O:9) LcrV is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.

**Purity:**

Greater than 80.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.