

Human AHA1 Recombinant Protein



RPPB5192

Product Information Protein Information

Product SKU:

RPPB5192

Accession:

O95433

Host:

Escherichia Coli.

Protein description:

AHA1 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 320 amino acids (19-337) and having a molecular mass of 36.1 kDa.

Appearance:

Sterile Filtered colorless solution.

Synonyms:

AHSA1, ACA1, p38, C14orf3, Activator of HSP90 ATPase-1.

Formulation:

The AHA1 protein solution contains 20mM Tris-HCl buffer (pH 8.0) and 2mM EDTA.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MATNVNNWHW TERDASNWST DKLKTLFLAV QVQNEEGKCE VTEVSKLDGE ASINNRKGKL IFFYEWSVKL
NWTGTSKSGV QYKGHVEIPN LSDENSDEV EISVSLAKDE PDTNLVALMKEEGVKLLREA MGIYISTLKT
EFTQGMILPT MNGESVDPVG QPALKTEERK AKPAPSKTQARPVGVKIPTC KITLKETFLT SPEELYRVFT
TQELVQAFTH APATLEADRG GKFHMDVGNV SGEFTDLVPE KHIVMKWRFK SWPEGHFATI TLTFIDKNGE
TELCMEGRGI PAPEEERTRQ GWQRYFEGI KQTFGYGARL