

Bovine TROVE2 Recombinant Protein



RPPB5221

Product Information Protein Information

Product SKU:

RPPB5221

Host:

Bovine tissues.

Protein description:

TROVE2 is a biotinylated natural protein having a molecular mass of 60kDa. TROVE2 is purified from bovine tissues by protein-chemical methods.

Appearance:

Sterile Filtered clear solution.

Synonyms:

60 kDa SS-A/Ro ribonucleoprotein, 60 kDa ribonucleoprotein Ro, RoRNP, 60 kDa Ro protein, Ro 60 kDa autoantigen, TROVE domain family member 2, Sjogren syndrome type A antigen, SS-A, Sjogren syndrome antigen A2, TROVE2, RO60, SSA2, RO-60.

Formulation:

TROVE2 is supplied in 20mM Tris/HCl buffer pH-7.5, 150mM NaCl, 2mM MgCl₂, 0.2mM DTT, 0.1mM PMSF and 20% glycerol.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.