

Phospholipase A2 Recombinant Protein



RPPB5264

Product Information Protein Information

Product SKU:

RPPB5264

Accession:

P00630

Host:

Escherichia Coli.

Protein description:

The E.Coli derived recombinant protein contains phospholipase P00630 bee venom protein epitopes, 26-162 amino acids. Phospholipase A2 P00630 Bee Venom is fused to a 6xHis tag at C-terminus

Synonyms:

Phospholipase A2, Phosphatidylcholine 2-acylhydrolase, Allergen Api m I, Api m 1.

Formulation:

20mM Tris-HCl, pH 7.2, 1.5M urea and 50% glycerol.

Purity:

Protein is >90% pure as determined by 10% SDS-PAGE.

Stability:

Phospholipase A2 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.