

Soybean P34 Recombinant Protein



RPPB5265

Product Information Protein Information

Product SKU:

RPPB5265

Host:

E.Coli.

Protein description:

The E.Coli derived recombinant protein contains Soybean P34 Protein full length, having an Mw of 30kDa. The protein is fused to a His tag at C-terminus.

Appearance:

Sterile filtered colorless solution.

Formulation:

50mM Tris-HCl, pH 8.0 and 60mM NaCl.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

Soybean P34 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles