

Soybean P34 GST Recombinant Protein



RPPB5266

Product Information Protein Information

Product SKU:

RPPB5266

Host:

Escherichia Coli.

Protein description:

The E.Coli derived GST Tag recombinant protein, 36kDa, contains Soybean P34 Protein epitopes 214-261, 351-379 amino acids.

Formulation:

50mM Tris-HCl, pH 8.0, 60mM NaCl, 10mM glutathione and 50% glycerol.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

Soybean P34 His although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.