

# Der F1 Recombinant Protein



RPPB5268

---

## Product Information    Protein Information

**Product SKU:**

RPPB5268

**Accession:**

P16311

**Host:**

E.Coli.

**Protein description:**

The E.Coli derived recombinant protein contains the Dermatophagoides farina Der F1 full length protein (a.a. 1-320) and fused to a 6 His Tag at C-terminus, having a total Mw of 36kDa, pI 5.88.

**Appearance:**

Sterile filtered colorless solution.

**Formulation:**

60mM NaCl, 50mM Tris-HCl pH 8.0 and 1.2M Urea.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

**Stability:**

Der F1 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.