

Aln G 4.0101 Recombinant Protein



RPPB5272

Product Information Protein Information

Product SKU:

RPPB5272

Accession:

O81701

Host:

Sf9 insect cells.

Protein description:

Recombinant Polcalcin Aln g 4 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 10,185 Dalton. Aln G 4.0101 is expressed with a 6xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Polcalcin Aln g 4, Calcium-binding pollen allergen Aln g 4, Aln g 4.

Formulation:

Aln G 4.0101 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.