

Bet v 2.0101 Recombinant Protein



RPPB5275

Product Information Protein Information

Product SKU:

RPPB5275

Accession:

P25816

Host:

Escherichia Coli.

Protein description:

Recombinant Profilin-1 produced in E.coli is a non-glycosylated, polypeptide chain having a calculated molecular mass of 15,625 Dalton. Bet v 2.0101 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Profilin-1, Allergen Bet v II, Pollen allergen Bet v 2, Bet v 2, BETVII, Bet v 2.0101.

Formulation:

Bet v 2.0101 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.