

Car b 1.0109 Recombinant Protein



RPPB5278

Product Information Protein Information

Product SKU:

RPPB5278

Accession:

B6RQR6

Host:

Escherichia Coli.

Protein description:

Recombinant Pollen allergen Car b 1 isoform produced in E.coli is a non-glycosylated, polypeptide chain having a calculated molecular mass of 18,757 Dalton. Car b 1.0109 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Pollen allergen Car b 1 isoform, Car b 1, Car b 1.0109.

Formulation:

Car b 1.0109 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.