

# Cor a 1.0103 Recombinant Protein



RPPB5279

---

## Product Information    Protein Information

**Product SKU:**

RPPB5279

**Accession:**

Q08407

**Host:**

Sf9 insect cells.

**Protein description:**

Recombinant Major pollen allergen Cor a 1 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 17,470 Dalton. Cor a 1.0103 purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Major pollen allergen Cor a 1 isoforms 5, 6, 11 and 16, Allergen Cor a I, Cor a 1, Cor a 1.0103.

**Formulation:**

Cor a 1.0103 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

**Purity:**

Greater than 80.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.