

Cor a 1.0401 Recombinant Protein



RPPB5280

Product Information Protein Information

Product SKU:

RPPB5280

Accession:

Q9SWR4

Host:

Sf9 insect cells.

Protein description:

Recombinant Cor a 1.0401 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 18,952 Dalton. Cor a 1.0401 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Major allergen Cor a 1.0401, CORA1.0401.

Formulation:

Cor a 1.0401 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.