

Phl p 1 Recombinant Protein



RPPB5281

Product Information Protein Information

Product SKU:

RPPB5281

Accession:

P43213

Host:

Sf9 insect cells.

Protein description:

Recombinant Pollen allergen Phl p 1 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 29,516 Dalton. Phl p 1 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Pollen allergen Phl p 1, Allergen Phl p I, Phl p 1, PHLPI.

Formulation:

Phl p 1 protein solution is supplied in 20mM HEPES pH-8.0, 200mM NaCl and 20% glycerol.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.