

Phl p 5 Recombinant Protein



RPPB5283

Product Information Protein Information

Product SKU:

RPPB5283

Accession:

O81342

Host:

Sf9 insect cells.

Protein description:

Recombinant Group V allergen Phl p 5.0203 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 27,582 Dalton. Phl p 5 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Group V allergen Phl p 5.0203, Phl p 5.

Formulation:

Phl p 5 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.