

Phl p 6 Recombinant Protein



RPPB5284

Product Information Protein Information

Product SKU:

RPPB5284

Accession:

P43215

Host:

Sf9 insect cells.

Protein description:

Recombinant Pollen allergen Phl p 6 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 13,162 Dalton. Phl p 6 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Pollen allergen Phl p 6, Allergen Phl p VI, Phl p 6, PHLPVI.

Formulation:

Phl p 6 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.