

Human UTI Recombinant Protein



RPPB5285

Product Information Protein Information

Product SKU:

RPPB5285

Protein description:

Ulinastatin is derived from human urine.

Host:

Human Urine.

Appearance:

Sterile Filtered White lyophilized (freeze-dried) powder.

Synonyms:

UTI, Ulinastatin, Urinary Trypsin Inhibitor.

Formulation:

Lyophilized from a (1mg/ml) solution containing no additives.

Solubility:

It is recommended to reconstitute the lyophilized UTI in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Stability:

Store the lyophilized UTI between 2-8°C, do not freeze. Upon reconstitution UTI should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

Biological Activity:

Human UTI has an activity of 2350IU/mg.