

# Phl p 12 Recombinant Protein



RPPB5286

---

## Product Information   Protein Information

**Product SKU:**

RPPB5286

**Accession:**

P35079

**Host:**

Sf9 insect cells.

**Protein description:**

Recombinant Phl p 12 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 15,607 Dalton. Phl p 12 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Profilin-1, Allergen Phl p 11, Pollen allergen Phl p 12, Phl p 12, PRO1, PHLPXI.

**Formulation:**

Phl p 12 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

**Purity:**

Greater than 80.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.