

# Fra e 1.0101 Recombinant Protein



RPPB5287

---

## Product Information    Protein Information

**Product SKU:**

RPPB5287

**Accession:**

Q7XAV4

**Host:**

Sf9 insect cells.

**Protein description:**

Recombinant Allergen Fra e 1.0101 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 17,794 Dalton. Fra e 1.0101 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Allergen Fra e 1.0101, Fra e 1.0101, Fra e 1.

**Formulation:**

Fra e 1.0101 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

**Purity:**

Greater than 80.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.