

Ara h 1.0101 Recombinant Protein



RPPB5288

Product Information Protein Information

Product SKU:

RPPB5288

Accession:

P43238

Host:

Sf9 insect cells.

Protein description:

Recombinant Allergen Ara h 1.0101 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 63,484 Dalton. Ara h 1.0101 is expressed with a 9xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Allergen Ara h 1, clone P41B, Allergen Ara h I, Allergen, Ara h 1.

Formulation:

Ara h 1.0101 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.