

Ara h 2.0201 Recombinant Protein



RPPB5289

Product Information Protein Information

Product SKU:

RPPB5289

Accession:

Q6PSU2

Host:

Sf9 insect cells.

Protein description:

Recombinant Allergen Ara h 2.0201 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 20,000 Dalton. Ara h 2.0201 is expressed with a 10xHis tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Conglutin-7, 2S protein 1, Seed storage protein SSP1, Seed storage protein SSP2, Allergen Ara h 2.

Formulation:

Ara h 2.0201 is supplied in 20mM HEPES buffer pH-7.9 and 6M Urea.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.