

MALD1 Recombinant Protein



RPPB5294

Product Information Protein Information

Product SKU:

RPPB5294

Accession:

Q9SYW3

Host:

Sf9 insect cells.

Protein description:

Recombinant MALD1 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 17,492 Dalton. MALD1 purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Major allergen Mal d 1, Ypr10 protein, MALD1, ypr10.

Formulation:

MALD1 is supplied in 20mM HEPES buffer pH-7.6, 250mM NaCl and 20% glycerol.

Purity:

Greater than 80.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.