

PfHSP70 Recombinant Protein



RPPB5296

Product Information Protein Information

Product SKU:

RPPB5296

Protein description:

The E.Coli derived recombinant Plasmodium falciparum contains the Plasmodium falciparum HSP70 protein epitopes 33-114 amino acids.

Appearance:

Sterile Filtered clear solution.

Formulation:

150mM Imidazole, pH 8.0, 150mM NaCl, 150mM Sodium Phosphate and 50% glycerol.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining) and RP-HPLC.

Stability:

PfHSP70 although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.