

Mumps Recombinant Protein



RPPB5300

Product Information Protein Information

Product SKU:

RPPB5300

Host:

Escherichia Coli.

Protein description:

Full length mumps viral nucleoprotein was expressed from E. coli, the recombinant protein was purified by affinity chromatography technique. The expressed recombinant mumps viral nucleoprotein contains 553 amino acids and a 6 x histidine tag was fused at its N-terminus, having a total Mw of 62kDa.

Appearance:

Sterile Filtered colorless solution.

Formulation:

Phosphate buffer, pH 7.4.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.