

Shiga Like Toxin 1 Recombinant Protein



RPPB5301

Product Information Protein Information

Product SKU:

RPPB5301

Host:

Escherichia Coli.

Protein description:

Recombinant Shiga Like Toxin-1 Subunit B from E.Coli O157:H7 produced in E.coli containing amino acids 2-90 of the Shiga Like Toxin-1 Subunit B. The Shiga Like Toxin 1 protein is fused to a 6xHis tag at its N-terminus and purified by proprietary chromatographic technique.

Formulation:

Phosphate buffered saline and 50mM arginine, pH 7.4.

Purity:

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.