

# Shiga Like Toxin 2 Recombinant Protein



RPPB5302

---

## Product Information    Protein Information

**Product SKU:**

RPPB5302

**Host:**

Escherichia Coli.

**Protein description:**

Recombinant Shiga Like Toxin-2 Subunit B from E.Coli O157:H7 produced in E.coli having an Mw of 22kDa. The Shiga Like Toxin 2 protein is fused to a 6xHis tag at its N-terminus and purified by proprietary chromatographic technique.

**Formulation:**

Shiga Like Toxin 2 protein solution contains PBS, 50mM arginine and 0.05% sodium azide.

**Purity:**

Protein is >95% pure as determined by 12% PAGE (coomassie staining).

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.