

SARS Mosaic S(M) Recombinant Protein



RPPB5304

Product Information Protein Information

Product SKU:

RPPB5304

Protein description:

The E.Coli derived 38 kDa recombinant mosaic protein contains the middle section of the Spike protein 408-470, 540-573 amino acids immunodominant regions.

Formulation:

25mM Tris-HCl, 0.4% sarcosyl, 0.25% Triton –100 and 50% glycerol.

Purity:

SARS Mosaic protein is >95% pure as determined by 10% PAGE (Coomassie staining).

Stability:

SARS Mosaic protein is shipped at ambient temperature. Upon arrival, store at -20C.