

SARS Mosaic S(C) Recombinant Protein



RPPB5305

Product Information Protein Information

Product SKU:

RPPB5305

Protein description:

The E.Coli derived 37 kDa recombinant mosaic protein contains the C-terminal t section of the Spike protein 1051-1076, 1121-1154, 1162-1190 amino acids immunodominant regions.

Formulation:

25mM Tris-HCl, 0.4% sarcosyl, 0.25% Triton –100 and 50% glycerol.

Purity:

SARS Mosaic Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

SARS Mosaic protein is shipped at ambient temperature. Upon arrival, store at -20C.