

SARS Core (340-390) Recombinant Protein



RPPB5306

Product Information Protein Information

Product SKU:

RPPB5306

Host:

Escherichia Coli.

Protein description:

The E.Coli derived 32kDa recombinant protein contains the Nucleocapsid core protein 340-390 amino acids immunodominant regions.

Appearance:

Sterile Filtered clear solution.

Formulation:

50mM Tris-HCl, 60mM NaCl, pH 8 and 50% glycerol.

Purity:

SARS Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

Protein is shipped at ambient temperature. Upon arrival, Store at -20°C.