

# SARS Nucleocapsid (1-49) Recombinant Protein



RPPB5307

---

## Product Information    Protein Information

**Product SKU:**

RPPB5307

**Protein description:**

The E.Coli derived 32 kDa recombinant protein contains the Nucleocapsid core protein 1-49 amino acids immunodominant regions. SARS Nucleocapsid protein 1-49 a.a. is fused to a GST tag.

**Formulation:**

50mM Tris-HCl, 60mM NaCl and 50% glycerol.

**Purity:**

SARS Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

SARS Core protein is shipped at ambient temperature. Upon arrival, store at -20C.