

# SARS SARS Core (1-49,192-220 a.a.) Recombinant Protein



RPPB5308

---

## Product Information    Protein Information

**Product SKU:**

RPPB5308

**Protein description:**

The E.Coli derived 35 kDa recombinant protein contains Nucleocapsid protein 1-49, 192-220 amino acids immunodominant regions. The 2 regions are separated with 3 glycine residues.

**Formulation:**

50mM Tris-HCl, pH-8, 60mM NaCl and 50% glycerol.

**Purity:**

SARS Core protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

SARS Core protein is shipped at ambient temperature. Upon arrival, store at -20C.