

SARS Nucleocapsid (2-422) Recombinant Protein



RPPB5310

Product Information Protein Information

Product SKU:

RPPB5310

Accession:

AAR87518.1

Host:

HEK293

Protein description:

The HEK293 derived recombinant protein contains the SARS Coronavirus Nucleoprotein, amino acids 2-422 fused to His tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

SARS Nucleocapsid protein solution is supplied in 0.2M NaCl & 50mM HEPES pH-8.

Purity:

Protein is >85% pure as determined SDS-PAGE.

Stability:

SARS nucleoprotein is shipped on ice packs. Upon arrival, Store at -20°C.