

# SARS Matrix Recombinant Protein



RPPB5311

---

## Product Information    Protein Information

**Product SKU:**

RPPB5311

**Protein description:**

The E.Coli derived 30kDa recombinant protein contains the Matrix protein 182-216 amino acids immunodominant regions.

**Formulation:**

50mM Tris-HCl, 60mM NaCl, and 50% glycerol.

**Purity:**

SARS Matrix Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

SARS Matrix protein is shipped at ambient temperature. Upon arrival, store at -20C.