

SARS Envelope Recombinant Protein



RPPB5312

Product Information Protein Information

Product SKU:

RPPB5312

Protein description:

The E.coli derived 34 kDa recombinant protein contains the N-terminus Enevelope protein 1-76 amino acids immunodominant regions.

Formulation:

50mM Tris-HCl, 60mM NaCl, pH 8 and 50% glycerol.

Purity:

SARS Envelope protein is >95% pure as determined by 10% PAGE(coomassie staining).

Stability:

SARS Envelope Protein is shipped at ambient temperature.Upon arrival, store at -20C.