

SARS MERS Recombinant Protein



RPPB5315

Product Information Protein Information

Product SKU:

RPPB5315

Accession:

AHX71946.1

Host:

HEK293

Protein description:

The HEK293 derived recombinant protein contains the SARS MERS Spike Glycoprotein S1, amino acids 18-725 fused to Fc tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

SARS MERS S1 protein solution is supplied in DPBS.

Purity:

Protein is >85% pure as determined SDS-PAGE.

Stability:

SARS MERS Protein is shipped on ice packs. Upon arrival, Store at -20°C.