

SARS Envelope Recombinant Protein



RPPB5317

Product Information Protein Information

Product SKU:

RPPB5317

Accession:

P0C6X7.1

Host:

Escherichia Coli.

Protein description:

The E.coli derived recombinant protein contains the Envelope protein 1-76 amino acids immunodominant regions fused to 6xHis tag at C-terminal.

Appearance:

Sterile Filtered clear solution.

Formulation:

SARS Envelope protein solution is supplied in PBS.

Purity:

Protein is >90% pure as determined SDS-PAGE.

Stability:

Protein is shipped on ice packs. Upon arrival, Store at -20°C.